Urine Processing Guide


Cup:



Used for collection of urine.

Urinalysis:

FILL BETWEEN THESE LINES



DO NOT REMOVE YELLOW CAP

TEST NAME/TEST CODE:

Urinalysis Complete (w/Microscopic)/**UCOMP**
Urinalysis Chemical/**UMAC**
Urinalysis Microscopic/**UMIC**

Urine Culture:

FILL TO THIS LINE



DO NOT REMOVE GRAY CAP

TEST NAME/TEST CODE:

Urine Culture (C&S, UC, UCX)/**CULURN**

REFLEX TESTING:

MUST SEND BOTH YELLOW AND GRAY TUP TUBE

Urinalysis (Reflex to Culture)/**UARFXC**
Urinalysis (Reflex to Microscopic)/**UARFXM**

CompuNet
CLINICAL LABORATORIES